

# Upright Laboratory Freezer (-40°C)

New stock coming – we need the room!

**Just 9,550**

**Usually \$14,626**

Code: SANYMDF-U443-PE

**PHCbi Model MDF-U443-PE Upright Biomedical-40°C Freezer** is ideally suited for the storage of plasma and DNA samples. The cascade refrigeration system ensures fast recovery from door openings. Direct fan forced cooling ensures superior chamber temperature uniformity.

### Key Features

- Alarms for hi/low temperature, power failure, filter, battery & fan motor replacement
- 5 polyethylene coated adjustable wire shelves (50kg load/shelf)
- Dry remote alarm contacts
- Door lock for security
- Stainless steel interior
- Access port; 40mm Ø
- Forced air circulation for superior chamber temperature uniformity
- Auto defrost has minimal effect on temperature of stored samples
- Manual defrost selectable.
- Two acrylic resin inner doors allow viewing of contents and maintain chamber temperature
- Easy access filter for simplified routine maintenance
- Stainless steel interior for durability & ease of cleaning
- 5 year warranty

**35% discount**

**Effective Capacity:** 426L

**Temperature Range:** -15°C to -40°C

**Exclusive Offer for a Limited Time. Contact Us to Order Today!**