

## An Exceptional, Affordable 406L Unit Suitable For A Wide Range Of Applications

**In Stock**

**Only \$9,100**

**Usually \$12,547**

Code: SANYMIR-554-PE

**Includes:**

5 adjustable height shelves



The wide variety of temperatures and lighting patterns that are essential in research, environmental studies and testing can now be accurately reproduced and controlled.

The stainless steel interior finish further enhances the durability of the unit whilst making it easy to clean.

- Wide temperature range:  
-10°C to +60°C with excellent uniformity
- Viewing window in outer door for sample observation
- Automatic and manual defrost
- Key lock for added security
- Light timer control (15W x 1) with optional light addition kit (MIR-L5) for adding 3 more fluorescent lights into chamber ceiling
- High and low-temperature alarms
- Access port for remote probes and power cords
- Independent over temperature protection
- 5-year warranty



Precise microprocessor temperature  
control with new LCD display

Ref. Code: SANY-MIR-554

Valid Until: 30 Apr 2025