

## 725L Upright -86°C Ultra-Low Temperature Freezer




**Only \$17,775**

**Usually \$25,545**

Code: SANYMDF-DU700ZH-PE

Tired of removing frost and ice from the inner door(s) of your ULT Freezer?

Then the FrostLess is the one for you, with improved insulation and a larger heating pipe frame, minimising frosting and maintenance.



- Air intake port for easy access to freezer chamber shortly after door closure
- Storage capacity:  
725 litres = 528 x 2" cryoboxes ( in optional racks)
- ECO "green gas" refrigerants and inverter compressors ensure excellent energy efficiency with reduced heat output
- Vacuum release port for faster access to freezer after door closure
- Inverter variable speed compressors use less energy by varying speed to match the required load
- Alarms for high and low temperature, door open and power failure
- USB port for downloading logged data
- Display keylock prevents adjustment of control panel by unauthorised personnel
- Compressor delay timer for staggered start up after power outage preventing circuit overload
- Castors for mobility and locating feet when in working position
- Keyed door locks for security
- Dry remote alarm contacts
- 5-year warranty

Ref. Code: SANY-Frostless

Valid Until: 30 Apr 2025