

# SMART

## Front Loading Autoclave



**Only \$50,910**

Usually \$56,571

**Model EH-200**

Code: PRIOEH200-SMART

With a cylindrical 200L chamber and  
2 non-tilting removable shelves

Chamber and door carry a  
10-year warranty



Biomaster  
Protected



Water  
Efficient



Energy  
Efficient



British  
Made

- Quickseal single-action door fitted with thermal and pressure locks to prevent opening at load temperature above 80 °C and pressures above 0.2 bar
- Operation up to 140C - 2.5 Bar
- 10 program memory
- Water cooled drain protection as standard
- Complying with AS/NZS standards
- Grade 316L stainless steel chamber
- TACTROL 3 microprocessor control system with simple and fully variable setting of process time, temperature and with graphic indication of cycle status
- Large cylindrical chamber within a compact frame
- Forced air cooling as standard helping to reduce cycle times

- Automatic timed & pulsed free steaming as standard for improved air removal and sterilising performance
- Steam generated within the chamber (electrically heated model)
- Complete with thermocouple entry and test ports
- Delayed start and media warming function as standard for ready to pour media at the start of the day
- Low water level sensor
- Epoxy coated panels and frame members are treated with an anti-bacterial agent which is effective against all bacteria and fungi including MRSA
- Balance of unit has a two year warranty
- Power Supply 400V 50Hz 10.5kW 15 amps

Ref. Code: PRIO-EH-200

Valid Until: 30 Apr 2025

