

Precision Incubation for Microtitre Plates & Petri Dishes

In Stock

Only \$5,525

Usually \$7,425

Code: TECPIM000510558

Includes a 3-year warranty

Experience precision, protection, and
ease of use with the Cole-Parmer®
Stuart™ INC-200D-M

The INC-200D-M is designed for uniform, reliable incubation, ensuring optimal sample conditions. With forced air circulation, preheated air flows evenly across the shelves, while stainless-steel water trays maintain humidity, preventing evaporation.

- **Even Heating** – Consistent temperature distribution for accurate results
- **Spacious & Accessible** – Stainless-steel chamber with 4 easy-slide shelves (fits 6-12 plates)
- **Built-in Safety** – Automatic temperature tracking & overheating protection
- **User-Friendly Display** – Bright LED screen for quick adjustments & alerts



Capacity: 20L

Temperature Range: Ambient +8°C to 80°C

Temperature Fluctuation (at 37°C): ± 0.5°C

Temperature Variation: ± 0.5°C

Ref. Code: COLE-200D

Valid Until: 31 May 2025