


No Wire Autoclavable Polypropylene Racks

Only \$45

Usually \$69

10mm - 13mm Tubes 90 Places	BELAF187460000	Red
	BELAF187450000	White
13mm - 16mm Tubes 60 Places	BELAF187470001	Blue
	BELAF187460001	Red
	BELAF187450001	White
16mm - 20mm Tubes 40 Places	BELAF187450002	White
25mm - 30mm Tubes 21 Places	BELAF187450004	White



A great alternative to wire racks, with price and space saving advantages over other plastic racks.

Polypropylene Test Tube Basket

Only \$69

Usually \$115

Code: BELAF187380000

Ideal for rinsing, transport and
storage of glassware.
Steam autoclavable at 121°C (250°F)
15.2cm x 15.2cm x 15.2cm



Ref. Code: Bel-Art_Sale

Valid Until: 31 Jul 2025

Clear Biohazard Disposal Bags

Only \$211

Pack of 100

Usually \$350

Code: BELAF131610009



Great economical solution for safe disposal of used pipettes, Petri dishes, multi-well cell culture plates, cell culture flasks, culture plates, membrane filters and other biohazard wastes.

Quaver Vibrating Spatula Blade

Only \$110

Usually \$168

Stainless Steel 19cm

Code: BELAH367140000



Make accurate additions of precious materials as needed. Rotating the toothed thumb wheel vibrates the curved tip spatula, giving you control over transfer from spatula to weighing dish.

Shatterproof Plastic Hydrometer

Only \$208

Usually \$320

Code: BELAB618000700

- Durable and shatterproof
- Easy to read black scale on yellow back helps avoid eye fatigue
- 255mm (10") long
- Lead free ballast
- Individually serialized



Combined Form Thermo-hydrometer

Only \$184

Usually \$282

Code: BELAB618210400

- Conforms to the accuracy requirements for ASTM thermo-hydrometers
- Used in industries such as food & beverage, chemical, pharmacology
- 381mm (15") long
- Lead free ballast and glass
- Temperature standardization: 60°F (15°C)



Ref. Code: Bel-Art_Sale

Valid Until: 31 Jul 2025