



ION-100A Ionizer

Eliminate Static. Maximize Accuracy.

Ideal for analytical work where every microgram matters.



Only \$2,295


Usually \$2,885

Code: OHAU30130303

Static charges can skew precision weighing especially with fine powders or small-volume samples.

The ION-100A neutralizes electrostatic interference, ensuring your OHAUS balance delivers consistently accurate results.

DC Corona Method



Key Features:

- **DC Corona Method:** Two electrodes work simultaneously to eliminate static between container and balance for accurate, stable readings
- **Fan-Free Operation:** Ionizes without disturbing delicate samples
- **Compact & Adjustable:** Fits neatly beside your balance with flexible positioning
- **Long-Lasting:** Durable, easily replaceable electrodes with a 15,000-hour lifespan

Ref. Code: OHAU-ION100A

Valid Until: 30 June 2025



0800 34 24 66
sales.nz@bio-strategy.dksh.com
www.bio-strategy.com
labshop.dksh.co.nz



Prices and discounts may be subject to change and only available while stocks last. Not valid with any other special offer or discounts and cannot be applied in retrospect. All prices in local currency excluding GST and local delivery charges. Products advertised may be for research use only (RUO) or in vitro diagnostic use (IVD). Please refer to manufacturer's product insert for further information. For product use outside of that specified on the manufacturer's labelling and/or instructions for use, Bio-Strategy does not take responsibility.