

## Easy CO<sub>2</sub> Incubator Decontamination with High Temperature Sterilisation



**Only \$10,245**

Usually \$13,661

5-year warranty

Code: SANYMCO-171AICUVDL-PE



Superior culturing and with the ease of reactive contamination control via high temperature sterilisation and pro-active contamination control via the inCu-saFe™ corrosion-resistant stainless steel copper enriched alloy interior!

### Key Features:

- Use of LED UV lamp reduces running costs.
- Full colour LCD touch panel provides clear parameter visibility.
- Dual wavelength single beam IR CO<sub>2</sub> sensor with continuous auto zero calibration.
- 4 adjustable height shelves.
- Direct heat and air jacket for environment stability.
- Programmable SafeCell UV system (isolated from cell cultures) irradiates airflow and humidity pan to prevent contamination without affecting cell cultures in-vitro.
- Humidity: 95%RH (±5%RH).
- Field reversible door.
- Peltier powered dew stick to draw condensation away from inner glass door.
- PID control system with colour LCD display.
- USB port for data transfer.
- Fully rounded corners and integrated shelf supports reduce contamination and simplify cleaning.

Ref. Code: SANY-CO<sub>2</sub>

Valid Until: 30 Jun 2025