

One of the most energy efficient ULT Freezers in it's class!



Only \$17,528

Usually \$25,041

729L

Capacity: 576
(2" cryoboxes)

7.7 kWh/day

Code: SANYMDF-DU702VH-PE



**30%
Discount**

** Pricing valid for current stock only*

The PHCbi VIP ECO Green Gas -86°C ultra-low temperature freezer uses natural refrigerants combined with inverter technology to provide efficient cooling and lower operating costs.

Key Features:

- Automatic vacuum relief port for quicker access to samples after door closure.
- Hot gas bypass combined with fivefold door gasket assists in minimising frosting and ice build up around door seal.
- Alarms for high and low temperature, door open and power failure.
- Three access ports.
- Three adjustable height stainless steel shelves creating four internal compartments.
- Use of natural (HC) refrigerants, inverter compressors and integrated electronics combine to lower the freezer running costs by up to 30% over conventional models.
- User friendly Large Colour Graphical Control Panel with digital display (can be used with gloved hand)
- USB port for downloading logged data.
- VIP Plus insulation: provides typically 24 % more storage capacity compared to conventionally insulated freezers with the same footprint.
- Castors for mobility and locating feet when in working position.
- Ezlatch door handle allows easy sample access and the keyed door lock provides security.
- Dry remote alarm contacts.
- Suitable for operation on 220/240V 50Hz 10 amp supply.
- Five Year Warranty.

Ref. Code: SANY-VIP-ECO

Valid Until: 30 Jun 2025

