

40L Front Loading Autoclave

Only \$24,760

Usually \$27,515

Code: PRIOB40-BASE

Includes:

316L stainless steel chamber
with 2 non-tipping shelves



Biomaster
Protected



Made in
Britain



Energy
Efficient



Water
Efficient

Big Performance in a Compact Unit!

Priorclave's benchtop autoclaves pack the performance and quality construction of a full-size, freestanding industrial autoclave into a footprint that fits any workspace. A quality research-grade benchtop autoclave offers all the flexibility, ease-of-use, and reliability that a small lab needs. Our compact tabletop units offer plug-n-play simplicity, but come with fully featured programmable control systems. They are perfect for media prep, liquid, and waste loads. A great fit for shared facilities and education.



Key Features:

- Quickseal single-action door fitted with thermal and pressure locks to prevent opening at load temperature above 80°C and pressures above 0.2 bar.
- Tactrol® 3 microprocessor control with simple and fully variable setting of process time and temperature with graphic indication of cycle status.
- Automatic timed free-steaming and pulsed free-steaming.
- Media Warming and Delayed Start functions keep autoclaved media at a usable temperature ready for immediate use.
- Low water level sensor.
- Two loading shelves.
- Epoxy coated panels and frame members are treated with an anti-bacterial agent which is effective against all bacteria and fungi including MRSA.
- Pressure vessel carries a 10 year warranty and is insurance approved.
- USB port
- Elements in chamber for ease of servicing & maintenance.
- Many options and accessories to select from.

Ref. Code: PRIO-B40

Valid Until: 30 Jun 2025

