



MICROFLEX® 93-260 Gloves

Maximum Protection, Minimal Compromise




Only \$240

Per Case

Case consists of
10 boxes / 50 gloves per box

Code	Size
ANSE93-260-S	Small
ANSE93-260-M	Medium
ANSE93-260-L	Large
ANSE93-260-XL	Xtra Large



PRODUCT HIGHLIGHT

1. Three-layer design:
Offering resistance to harsh chemicals
2. Thinnest chemical-resistant disposable glove:
For tactility and dexterity
3. Low AQL: 0.65 AQL for allowable pinholes means
exceptional barrier integrity

Get in touch with our team for chemical
compatibility, samples & bulk pricing

Engineered for demanding environments,
the Ansell 93-260 offers exceptional
chemical resistance in one of the thinnest
designs on the market.

Its three-layer nitrile/neoprene barrier
protects against acids, bases, and a broad
range of solvents, while maintaining
superior tactile sensitivity and dexterity.

The extended cuff adds extra forearm
protection, and the textured grip ensures
control in both wet and dry conditions.

Ideal for laboratory, pharmaceutical, and
industrial applications, the 93-260 delivers
comfort, durability, and confidence —
even during prolonged wear.

Ref. Code: ANSE-93-260

Valid Until: 31 Aug 2025



0800 34 24 66
sales.nz@bio-strategy.dksh.com
www.bio-strategy.com
labshop.dksh.co.nz



Prices and discounts may be subject to change and only available while stocks last. Not valid with any other special offer or discounts and cannot be applied in retrospect. All prices in local currency excluding GST and local delivery charges. Products advertised may be for research use only (RUO) or in vitro diagnostic use (IVD). Please refer to manufacturer's product insert for further information. For product use outside of that specified on the manufacturer's labelling and/or instructions for use, Bio-Strategy does not take responsibility.