

**ESCO**  
WORLD CLASS. WORLDWIDE.

**NordicSafe®**



**4 In Stock**

**Class II 1.2m  
Microbiological Safety Cabinet  
with motorised sash**

**Only \$14,800**

**Usually \$19,592**

Code: ESCONC2-4L8

(Cabinet only)

Includes UV light & 2 power sockets as standard

Optional height adjustable stand

**\$995** (usually \$1,515)

Code ESCOSTC-4CO

Simply press the Up arrow to raise the sash to the correct working height (200 mm) and activate the blower fan and light.

- Ergonomic sloping front combined with tempered glass side panels allows entry of more natural light
- Airflow sensor monitors real-time airflow for safety and alerts the user if airflow is insufficient
- Motorised sash ensures sash travels to working height automatically and automatically activates internal lighting and airflow
- Night set back mode reduces power consumption by a further 60%
- ULPA downflow and exhaust filters with 99.999% efficiency, 10X better than HEPA filter (99.99%) improves safety and delivers ISO Class 3 air cleanliness
- Sentinel Gold microprocessor controller with digital display shows cabinet status
- Audible and visual alarms
- Raised armrest external to cabinet for comfortable working position
- Multi-piece stainless steel worksurface is easy to remove for cleaning
- Certified by PHE (Public Health England) for conformance to EN 12469
- Airflow sensor monitors real-time airflow for safety and alerts the user if airflow is insufficient
- Removable easy-to-clean paper catch
- Energy efficient DC ECM Motor (with 70% energy savings compared to AC motor)
- ISOCIDE™ Powder Coat – painted surfaces have silver ion impregnated coat to inhibit microbial growth
- 4-year warranty

Ref. Code: ESCO-Nordic

Valid Until: 31 Jul 2025

 **bio-strategy**  
Part of DKSH Group



0800 34 24 66  
sales.nz@bio-strategy.dksh.com  
www.bio-strategy.com  
shop.bio-strategy.com

 **DKSH**

Prices and discounts may be subject to change and only available while stocks last. Not valid with any other special offer or discounts and cannot be applied in retrospect. All prices in local currency excluding GST and local delivery charges. Products advertised may be for research use only (RUO) or in vitro diagnostic use (IVD). Please refer to manufacturer's product insert for further information. For product use outside of that specified on the manufacturer's labelling and/or instructions for use, Bio-Strategy does not take responsibility.