

## UP200St with VialTweeter for Simultaneous Sample Preparation

**Only \$14,830**

**Usually \$16,475**

Code: VialTweeter-Set

**Includes:**

- 1 x Ultrasonic processor
- 1 x VialTweeter Sonotrode



The VialTweeter, made from titanium, is a lab homogenizer that allows simultaneous sonication of up to 10 vials (Ø11mm) without cross-contamination or sample loss, making it ideal for preparing multiple small volume samples.

It uses high power ultrasound to achieve processes such as homogenization, emulsification, dispersing, and degassing by introducing ultrasonic cavitation through the vessel's wall

Combined with the digital ultrasonic processor UP200St, it ensures even and reproducible sonication results across all samples.

### Advantages at a glance

- Intense sonication of up to 10 vials simultaneously
- Indirect sonication at high ultrasonic intensity through vessel wall into the sample. Indirect sonication avoids cross-contamination and sample loss
- Reproducible results due to adjustable and controllable sonication amplitude
- The VialPress enables you to sonicate larger tubes
- Adjustable pulse mode from 0 to 100%
- Autoclavable



Ref. Code: HIEL-UP200ST

Valid Until: 31 Aug 2025