



Only \$3,400

Usually \$4,290

Code: IKAA2900005

Batch mill for 2 different grinding procedures:

Impact grinding of hard, brittle or non-elastic grinding materials with high-grade stainless steel beater.

This beater can be used for a Mohs hardness up to 6 (included).

Cutting grinding for pulverizing soft, fibrous materials with a cutting blade (not included).

- Moist and gluey materials can be pulverized by adding water
- Grinding chamber made of Tefcel (ETFE, glass fiber-reinforced) with stainless steel inlet (1.4571), useful volume 80 ml (incl. with delivery). For embrittlement of grinding materials with liquid nitrogen in the grinding chamber
- Optionally, a 250 ml grinding chamber is available

Optional accessories available.
Contact us for more information and pricing

Ref. Code: IKAA-A11

Valid Until: 31 Aug 2025