

Airstream®



ESCO
WORLD CLASS. WORLDWIDE.

Gen3 Horizontal Laminar Flow Cabinet



Only \$6,640

Usually \$8,895

Code: ESCOLHG-4AG-F8

Includes: UV light, front cover &
2 power outlets

Why is the Airstream 1.2m Horizontal Laminar Flow Cabinet so popular?

It is one of the quietest and safest cabinet choices for sample protection and meets many international standards.

- ULPA filter (10x filtration efficiency than of HEPA filter) delivers ISO Class 3 air cleanliness
- Noise level: only 55.8dB(A) ; Interior lighting: 1,394 Lux
- When not in use: Standby mode provides more energy savings and maintains an ISO Class 5 work zone
- Stainless steel spill retaining worksurface for ease of cleaning
- ISOCIDE™ Powder Coat – painted surfaces have silver ion impregnated coat to inhibit microbial growth
- 3-year warranty

Optional adjustable height stand

Code: ESCOSTC-4AO

Only \$995 | Usually \$1,515

Applications:

- Plant Tissue Culture
- Media Plate Preparation
- Electronics Inspection
- Medical Device Assembly

Ref. Code: ESCO-LHG

Valid Until: 31 Aug 2025

 **bio-strategy**
Part of DKSH Group



0800 34 24 66
sales.nz@bio-strategy.dksh.com
www.bio-strategy.com
labshop.dksh.co.nz



Prices and discounts may be subject to change and only available while stocks last. Not valid with any other special offer or discounts and cannot be applied in retrospect. All prices in local currency excluding GST and local delivery charges. Products advertised may be for research use only (RUO) or in vitro diagnostic use (IVD). Please refer to manufacturer's product insert for further information. For product use outside of that specified on the manufacturer's labelling and/or instructions for use, Bio-Strategy does not take responsibility.