

38m x 10cm

**Only \$65**

Usually \$90

Code: MARI7403810

75m x 5cm

**Only \$68**

Usually \$95

Code: MARI7407505



From sealing flasks and petri dishes to preventing evaporation and contamination, Parafilm® has become the go-to solution for everyday lab tasks. Its flexibility (clings around irregular shapes and surfaces), stretchability (up to 200%), and unbeatable seal make it a staple on every benchtop—because no lab feels complete without a roll within arm's reach.

### Storage:

Parafilm® M can be stored for at least 3 years at a temperature between 7 °C and 32°C and relative humidity of 50 %.

Ref. Code: MARI-PARA

Valid Until: 31 Aug 2025