



Endeavor™ 5000 Microplate Shaker

Only \$1,750

Usually \$2,265

Code: OHAU30770923



Motion: Orbital, 3 mm

Capacity: 4 Microplates or 2 Micro-tube Racks

Speed Range: 100 rpm – 1,200 rpm

Warranty: 2 years

Construction

Durable chemical-resistant polymer housing, sturdy easy-to-clean aluminum standard tray with attachment points for accessories, IP21 Ingress Protection, robust triple eccentric drive for stable shaking, brushless DC motor, control knobs, on/off switch

Design Features

Visual and audio alarm when timer counts down to zero or an error occurs, small footprint to save bench space, advanced microprocessor control slowly ramps the shaker tray up to the desired speed and avoids splashing, easy to attach optional accessory platform


Optional Accessories



Micro Tube Rack

Holds 70 x 1.5mL to 2mL microtubes

Code: OHAU30400114



Deep Plate Kit

To fix 4x deep well plates to shaker

Code: OHAU30772773

Ref. Code: OHAU-Endeavor

Valid Until: 31 Aug 2025



Part of DKSH Group



0800 34 24 66
sales.nz@bio-strategy.dksh.com
www.bio-strategy.com
labshop.dksh.co.nz



Prices and discounts may be subject to change and only available while stocks last. Not valid with any other special offer or discounts and cannot be applied in retrospect. All prices in local currency excluding GST and local delivery charges. Products advertised may be for research use only (RUO) or in vitro diagnostic use (IVD). Please refer to manufacturer's product insert for further information. For product use outside of that specified on the manufacturer's labelling and/or instructions for use, Bio-Strategy does not take responsibility.