

238L Cooled Incubator

Only \$6,344

Usually \$7,931

Supplied with 5 shelves

Code: SANYMIR-254-PE

In Stock

The affordable solution to your
cooled incubation needs!

An exceptional unit suitable for a wide range of
applications requiring a -10°C to +60°C environment.



- Key lock for added security
- Viewing window in outer door for sample observation
- Light timer control (15W x 1) with optional light addition kit (MIR-L5) for adding 3 more fluorescent lights into chamber ceiling
- Glass viewing window in door for easy viewing of product
- Access port for remote probes and power cords
- Precise microprocessor temperature control with new LCD display
- Automatic and manual defrosting
- Independent over temperature protection
- 5-year warranty

Ref. Code: SANY-MIR254

Valid Until: 31 Aug 2025