




Plant Growth and Environmental Chamber



Only \$14,702

Usually \$19,603

Code: SANYMLR-352-PE

Only \$18,820

Usually \$25,094

Code: SANYMLR-352H-PE

Supplied with 5 shelves, access port &
a 5-year warranty

The MLR-352-PE allows control of temperature and lighting, while
the MLR-352H can also control humidity.

- **Temperature range:**
0 to 50°C with lights off. +5 to 50°C with lights on
- **Lighting:**
to 20,000 lux via 15 x 40 watt fluorescent tubes in 6 increments
- **Humidity:**
60% to 90% RH lights off/ 55% to 85% RH lights on
- Suitable for plant, tissue and insect culturing
(insect culturing will require frequent maintenance)
- 10 program storage and each program can have up to 12-steps
with repeat 1-98 times or unlimited
- Alarms for temperature (upper & lower), sensor, filter and door open
- Light banks in front and side panels
- Microprocessor PID and refrigeration capacity control with graphic
LCD display
- Data Logging function: up to 14 days in 6 minute intervals
- Automatic and manual defrost

Ref. Code: SANY-Growth

Valid Until: 31 Aug 2025