

Only \$8,795

Usually \$14,696

Capacity: 426L with
5 adjustable shelves

Temp Range: -15°C to -40°C

Code: SANYMDF-U443-PE



Ideally suited for the storage of plasma and DNA samples.

The cascade refrigeration system ensures fast recovery from door openings.

Direct fan forced cooling ensures superior chamber temperature uniformity.

- Stainless steel interior for durability & ease of cleaning
- Alarms for hi/low temperature, power failure, filter, battery & fan motor replacement
- Dry remote alarm contacts
- Door lock for security
- Auto defrost has minimal effect on temperature of stored samples
- Manual defrost selectable
- Two acrylic resin inner doors allow viewing of contents and maintain chamber temperature
- Easy access filter for simplified routine maintenance
- 5-year warranty

Ref. Code: SANY-Plasma

Valid Until: 31 Aug 2025