

## Front Loading 320L Autoclave

### EH320-BASE

suitable for non-vacuum  
sterilisation runs  
5 program memory

**Only \$59,385** | Usually \$65,984

Code: PRIOEH320-BASE

### EH320-SMART

vacuum options: pre-cycle, cooling &  
drying but can be used without  
10 program memory

**Only \$67,000** | Usually \$74,446

Code: PRIOEH320-SMART

Export Crate

Code PRIOPC/AEC/020 **Only \$2,957** | Usually \$3,286



Many options and accessories available  
including loading trolleys, loading cassettes,  
printer, wandering probe, discard containers,  
water softener etc.

Please contact us for a site visit to  
assist in configuring your autoclave  
to your requirements.

The grade 316L stainless steel chamber is durable and the  
unit is inherently safe with thermal and pressure interlocks.

- Operation up to 140°C / 2.5 Bar
- Chamber Dimensions: 630mm Ø x 1,025mm deep
- Water cooled drain protection as standard
- Quickseal single-action door fitted with thermal and pressure locks to prevent opening at load temperature above 80°C and pressures above 0.2 bar
- TACTROL 3 microprocessor control system with simple and fully variable setting of process time, temperature and with graphic indication of cycle status
- USB port for logging, remote diagnosis and updating
- Automatic water filling
- Low water level sensor
- Forced air cooling as standard helping to reduce cycle times
- Automatic timed and pulsed free steaming as standard for improved air removal and sterilising performance
- Steam generated within the chamber (electrically heated model)
- Complete with thermocouple entry and test ports
- Delayed start and media warming function as standard for ready to pour media at the start of the day
- Epoxy coated panels and frame members are treated with an anti-bacterial agent which is effective against all bacteria and fungi including MRSA
- Pressure vessel carries a 10-year warranty and is insurance approved

Ref. Code: PRIO-EH320

Valid Until: 31 Aug 2025