

MaxyClear Snaplock Non-Sterile Microtubes



AXYGMCT-200-C

Size 0.6ml

Only \$37

Usually \$49

Code: AXYGMCT-060-C
[1,000 items]

Size 1.5ml

Only \$20

Usually \$23

Code: AXYGMCT-150-C
[500 items]

Size 1.75ml

Only \$18

Usually \$21


Code: AXYGMCT-175-C
[500 items]

Size 2.0ml

Only \$21

Usually \$26

Code: AXYGMCT-200-C
[500 items]



AXYGMCT-150-C

- Will withstand up to 16,700 RCF
- Snaplock cap assures an efficient seal
- Thin membrane in the cap centre allows easy access to the sample by syringe needle
- RNase-/DNase-free and nonpyrogenic for total sample integrity
- Frosted cap surfaces and side panels for labelling and writing

Axygen Microtubes - Precision You Can Trust

Engineered for reliability and performance, Axigen Microtubes offer superior clarity, secure sealing, and maximum sample integrity.

With options for sterile, RNase-/DNase-free, and low-retention formats, these tubes are ideal for critical molecular biology, PCR, and storage applications.

Designed to withstand rigorous centrifugation and freeze-thaw cycles, Axigen delivers the quality and consistency trusted by labs worldwide.

Ref. Code: AXYG-MCT

Valid Until: 30 Sep 2025



0800 34 24 66
sales.nz@bio-strategy.dksh.com
www.bio-strategy.com
labshop.dksh.co.nz



Prices and discounts may be subject to change and only available while stocks last. Not valid with any other special offer or discounts and cannot be applied in retrospect. All prices in local currency excluding GST and local delivery charges. Products advertised may be for research use only (RUO) or in vitro diagnostic use (IVD). Please refer to manufacturer's product insert for further information. For product use outside of that specified on the manufacturer's labelling and/or instructions for use, Bio-Strategy does not take responsibility.