

# HM-200 Electromantle



**Please Note:** the Electrothermal Electromantle Models have been rebranded Cole-Parmer under the HM-200 range. All the same great features and specifications remain, they just have a new and improved look.

Designed with flexibility and ease of use in mind, the HM-200 series heating mantle is a reliable workhorse for your heating applications in the laboratory.

Manufactured to the highest quality, each mantle is made with extreme care and is hand sewn to assure its reliability, quality and safety.

Key safety features include: built-in energy regulators, indicator lamps, double fuses and excellent insulation.

- Use with round bottom flasks from 50mL to 5L
- Maximum element temperature 450°C
- Chemically resistant polypropylene outer casing
- Supplied with bracket for 1.3cm diameter support rods (support rods not included)
- 3-year Warranty

Code	Capacity	Usually	Now
36002-02	50mL	\$1,075	<b>\$940</b>
36002-06	100mL	\$1,150	<b>\$1,005</b>
36002-10	250mL	\$1,160	<b>\$1,015</b>
36002-14	500mL	\$1,370	<b>\$1,195</b>
36002-18	1,000mL	\$1,360	<b>\$1,190</b>
36002-22	2,000mL	\$2,370	<b>\$2,070</b>
36002-26	3,000mL	\$2,245	<b>\$1,960</b>
36002-30	5,000mL	\$2,135	<b>\$1,865</b>

Ref. Code: COLE-HM-200

Valid Until: 30 Sep 2025