

# Class B Blue Burette

Reliable, easy to use, and now more affordable than ever

**Only \$45**

**Usually \$80**

Code: BIBBBRB266P



The DWK MBL® Class B 50 mL burette with modular PTFE stopcock is a practical choice for everyday lab work.

Whether you're teaching basic titrations, preparing solutions, or running routine tests, this burette offers all the functionality you need without the high price tag.

Thanks to current stock clearance pricing, it's a great time to stock up on trusted DWK quality at exceptional value.

- Complies with ISO 385
- Burette manufactured from borosilicate glass 3.3
- With fluted funnel top to make filling easier and safer
- PTFE stopcock for grease-free performance
- Modular design with detachable burette tube, stopcock, and jet for easy replacement
- Easy to read blue enamel graduations

Ref. Code: DWK-B266P

Valid Until: 30 Sep 2025