

# Zefon Cassette

Test for asbestos presence in  
the air with PCM cassettes



**Only \$93**

**Usually \$115**

25mm PCM Cassette with  
Green Grid Filters (50pk)  
Code: Z008GGA

Each 25-mm, three-piece ABS plastic cassette housing is injection molded in-house and is individually assembled and preloaded under stringent quality control standards, ensuring you get accurate results.

Once assembled, each cassette undergoes an inspection for the correct filter/pad combination before being shipped out.

Random lot sample testing is performed daily by an independent laboratory to ensure low background fiber counts.

Meets NIOSH, OSHA, and AHERA standards for asbestos air monitoring.

- PCM cassettes come preloaded, ready for use
- Offers fast turnaround and low cost
- Each cassette marked with a unique serial number
- 2" static conductive extension cowl
- Mixed cellulose ester membrane
- 0.8 µm pore size

Ref. Code: ENEX-Zefon

Valid Until: 30 Sep 2025