



# Vortex Mixer

Robust, reliable, and great value!

**Only \$425**

**Usually \$685**

Code: OHAU30392118

Each unit includes a cup head and a 3-inch accessory platform, accommodating a wide range of tubes and vessels.

The Ohaus Analogue Mini Vortex Mixer features variable speed control, enabling everything from gentle shaking to high-speed, vigorous vortexing of samples.

Choose between convenient touch mode for quick mixing or continuous mode for hands-free operation with a range of accessories. Its rugged, stable design keeps the unit firmly in place, even during vigorous use.

- Control: Analogue
- Orbit: 4.9 mm
- Speed Range:  
300 rpm – 2800 rpm
- 2-years Warranty



A wide range of option accessories also available including holders, microplates, foam inserts and more

Give us a call to discuss your mixing needs



Ref. Code: OHAU-Vortex

Valid Until: 31 Aug 2025



Part of DKSH Group



0800 34 24 66  
[sales.nz@bio-strategy.dksh.com](mailto:sales.nz@bio-strategy.dksh.com)  
[www.bio-strategy.com](http://www.bio-strategy.com)  
[labshop.dksh.co.nz](http://labshop.dksh.co.nz)



Prices and discounts may be subject to change and only available while stocks last. Not valid with any other special offer or discounts and cannot be applied in retrospect. All prices in local currency excluding GST and local delivery charges. Products advertised may be for research use only (RUO) or in vitro diagnostic use (IVD). Please refer to manufacturer's product insert for further information. For product use outside of that specified on the manufacturer's labelling and/or instructions for use, Bio-Strategy does not take responsibility.