



Compact Affordable -30°C

Biomedical Chest Freezers



Only \$2,841

Usually \$4,059

Code: SANYMDF-137-PE

138L including
2 baskets



**Limited
Stock Available**

Only \$3,194

Usually \$4,563

Code: SANYMDF-237-PE

221L including
3 baskets



- Temperature range: -20°C to -30°C
- Hi/Lo temperature & power failure alarms & dry remote alarm contacts
- Hermetic compressor with low energy consumption
- User friendly controller with digital display
- Display keylock prevents adjustment of control panel by unauthorised personnel

- Compressor delay timer for staggered start up after power outage preventing circuit overload
- Castors for mobility & locating feet for working position
- Quiet in operation : 41dB (A)
- Keyed lid lock for security
- Suitable for operation on 220/240V 50Hz 10 amp supply
- 5 year warranty

Ref. Code: SANY-Chest

Valid Until: 31 Aug 2025



Part of DKSH Group



0800 34 24 66
sales.nz@bio-strategy.dksh.com
www.bio-strategy.com
labshop.dksh.co.nz



Prices and discounts may be subject to change and only available while stocks last. Not valid with any other special offer or discounts and cannot be applied in retrospect. All prices in local currency excluding GST and local delivery charges. Products advertised may be for research use only (RUO) or in vitro diagnostic use (IVD). Please refer to manufacturer's product insert for further information. For product use outside of that specified on the manufacturer's labelling and/or instructions for use, Bio-Strategy does not take responsibility.