



PHCbi



TwinGuard 360L -86°C Freezer

For those who appreciate the
value of their samples!

Only \$13,995

Usually \$21,911

Code: SANYMFD-DU302VX-PE

Pricing valid for current stock only

Capacity of 360L
(240 x 2" cryoboxes)

Offering enhanced safety due to the dual cooling system. Two independent refrigeration systems provide the best redundancy you can have! – if one refrigeration circuit fails – the other circuit will maintain -70°C.

- User friendly Colour Graphical Control Panel with digital display
- USB port for downloading logged data
- Filter-less design reduces the time for routine maintenance
- Quiet in operation – only 52dB(A)
- Hot gas bypass combined with fivefold door gasket minimises frosting and ice build up around door seal
- Automatic vacuum relief port for quicker access to samples after door closure
- Alarms for high and low temperature, door open and power failure
- Intelligent ECO Mode: optimising compressor operation to match the freezer's running condition & safe electricity costs
- VIP Plus insulation: provides typically 24% more storage capacity compared to conventionally insulated freezers with the same footprint
- Castors for mobility and locating feet when in position
- Keyed door locks for security
- Dry remote alarm contacts
- 5-year warranty

Ref. Code: SANY-TWIN302

Valid Until: 30 Sep 2025



0800 34 24 66
sales.nz@bio-strategy.dksh.com
www.bio-strategy.com
labshop.dksh.co.nz



Prices and discounts may be subject to change and only available while stocks last. Not valid with any other special offer or discounts and cannot be applied in retrospect. All prices in local currency excluding GST and local delivery charges. Products advertised may be for research use only (RUO) or in vitro diagnostic use (IVD). Please refer to manufacturer's product insert for further information. For product use outside of that specified on the manufacturer's labelling and/or instructions for use, Bio-Strategy does not take responsibility.