



PHCbi



Frostless 725L -86°C Freezer

Only \$19,150

Usually \$25,545

Includes 3 shelves & 5-year warranty

Code: SANYDU700ZH-PE

Pricing valid for current stock only

**Imagine a ULT Green Gas “ECO” Freezer
with reduced maintenance.
Then the Frostless is your answer!**

- Storage capacity: 725L = 528 x 2” cryoboxes (in optional racks)
- Quiet in operation – only 52dB(A)
- VIP Plus insulation: provides typically 24% more storage capacity compared to conventionally insulated freezers with the same footprint
- Air intake port for easy access to freezer chamber shortly after door closure
- Inverter variable speed compressors use less energy by varying speed to match the required load
- Alarms for high and low temperature, door open and power failure
- User friendly Colour Graphical Control Panel with digital display
- USB port for downloading logged data
- Display keylock prevents adjustment of control panel by unauthorised personnel
- Compressor delay timer for staggered start up after power outage preventing circuit overload
- Castors for mobility and locating feet when in position
- Keyed door locks for security
- Dry remote alarm contacts

Denser packing/insulation in the inner doors combined with a larger hot gas bypass near the door gasket reduces frosting & ice build-up

ECO “green gas” refrigerants and inverter compressors ensure excellent energy efficiency with reduced heat output

Ref. Code: SANY-Frostless725

Valid Until: 30 Sep 2025

 **bio-strategy**
Part of DKSH Group



0800 34 24 66
sales.nz@bio-strategy.dksh.com
www.bio-strategy.com
labshop.dksh.co.nz

 **DKSH**

Prices and discounts may be subject to change and only available while stocks last. Not valid with any other special offer or discounts and cannot be applied in retrospect. All prices in local currency excluding GST and local delivery charges. Products advertised may be for research use only (RUO) or in vitro diagnostic use (IVD). Please refer to manufacturer's product insert for further information. For product use outside of that specified on the manufacturer's labelling and/or instructions for use, Bio-Strategy does not take responsibility.