

Looking to improve the efficiency of your formulation & synthetic processes?

We've got exciting news for you – Büchiglas, a powerful player in chemical reactor technology is now exclusively distributed by DKSH.

Here's what's coming your way:

- High-quality metal or glass reactor systems for chemical synthesis
- Available in lab scale configurations up to production plants
- Customisable configurations to suit every customer requirement
- Family owned, globally trusted, focussing on sustainability & continuous improvement



miniPilot: glass reactor for chemical synthesis with configurations ranging from 5, 10 and 15L.
 Max operating pressure of 1bar and up to 200C



Polyclave: small pressure reactor capable of working with volumes between 0.25 – 5.0L.
 Pressure rated up to 300bar and 250C

Interested to learn more?
 Our team are here ready to help

Ref. Code: BGLA-Reactor

Ref: Nov 2025