


STOCK CLEARANCE

Designed For Laboratory Culturing Applications



Fan Forced, Ambient + 7.5°C to 100°C
Supplied with 2 shelves



Pre-Heat Chamber Technology

4-zone heated air jacket for maximum
temperature uniformity

110 Litre

Only \$3,400

Usually \$5,210

Code: ESCOIFA-110-8

240 Litre

Only \$4,985

Usually \$8,463

Code: ESCOIFA-240-8

170 Litre

Only \$3,995

Usually \$6,847

Code: ESCOIFA-170-8

Order Today while stocks last!

Find more great
ESCO deals here →



- Isocide™ Antimicrobial Powder Coating: Eliminates 99.9% of surface bacteria within 24 hours of exposure
- Single-piece stainless steel chamber with directly mounted shelves making incubator easy to clean
- Permanently lubricated maintenance-free fan
- Air flow slider to adjust exchange rate of air
- Microprocessor PID control with 10 program memory

Ref. Code: ESCO-IFA

Valid Until: 31 Jan 2026