

# microman<sup>®</sup>

E

## MICROMAN E Bundle

 **GILSON<sup>®</sup>**

**\$555**

**\$999**

1 x MICROMAN E valued @ \$513  
1 x Capillary Piston valued @ \$143

1 x MICROMAN E valued @ \$513  
1 x Capillary Piston valued @ \$663



Contact your local  
Sales Rep for full details



Watch MICROMAN E in action

### MICROMAN E

FD10001	1-10ul
FD10002	3-25ul
FD10003	20-50ul
FD10004	10-100ul
FD10005	50-250ul
FD10006	100-1000ul

### Capillary Piston

F148412	CP10	192 items	\$143
F148560	CP1000	182 items	
F148012	CP25	960 items	\$663
F148013	CP50		
F148314	CP100		
F148014	CP250		

### ERGONOMIC PIPETTING FOR NON-AQUEOUS LIQUIDS

MICROMAN E positive-displacement pipettes, along with capillary pistons, offer the highest precision when ergonomically pipetting non-aqueous samples, such as viscous or volatile liquids. Capillary pistons eliminate the air cushion between the sample and disposable piston, ensuring the accuracy of the pipette is not affected by the physical properties of the sample.

#### Triple Protection

Positive-displacement technology is the ideal solution to prevent from aerosols contamination or sample carry-over, protecting the pipette, the user, and the sample.

#### Capillary Pistons

##### 100% Protection Against Contamination

The disposable piston is located inside the capillary tip so users can be assured that no harmful samples will contaminate the pipette.

- No risk of chemical contamination
- Slender profile—fits most tubes

Valid Until: 31 Mar 2026

Ref. Code: GILS-Bundle

0800 951 010 | customerservicetecnz@dksh.com | labshop.dksh.co.nz



Prices and discounts may be subject to change and only available while stocks last. Not valid with any other special offer or discounts and cannot be applied in retrospect. All prices in local currency excluding GST and local delivery charges. Products advertised may be for research use only (RUO) or in vitro diagnostic use (IVD). Please refer to manufacturer's product insert for further information. For product use outside of that specified on the manufacturer's labelling and/or instructions for use, Bio-Strategy does not take responsibility.

Delivering Growth

 **bio-strategy**  
Part of DKSH Group

 **DKSH**