

Biomedical Freezer 690L

Preservation
sub-zero



Only \$7,406

Usually \$10,581

Code: MDF-U731M-PE

**2 Units Only
at this great price**

Order Yours Today!



The MDF-U731M is a large-capacity laboratory freezer with high freezing power.

Specifically designed to offer a larger capacity freezer with the uniformity, stability, and reliability required in the life science market, chosen by researchers for its unrivalled performance and overall design.

- Selectable temperature range from -20°C to -30°C
- Very quiet in operation – only 42dB(A)
- Direct cooling system with manual push-button defrost
- Supplied with 4 adjustable polyethylene-coated wire shelves
- External Dimensions: 770 x 830 x 1,955mm W x D x H
- Alarms for high and low temperature, power failure, door open, and compressor temperature
- Two access ports.
- Castors for mobility and locating feet when in working position
- Door lock with padlock compatibility
- Dry remote alarm contacts
- Suitable for operation on 220/240V 50Hz 10 amp supply
- 5-year warranty

Valid Until: 28 Feb 2026

Ref. Code: SANY-U731M

0800 34 24 66 | sales.nz@bio-strategy.dksh.com | www.bio-strategy.com | labshop.dksh.co.nz

Prices and discounts may be subject to change and only available while stocks last. Not valid with any other special offer or discounts and cannot be applied in retrospect. All prices in local currency excluding GST and local delivery charges. Products advertised may be for research use only (RUO) or in vitro diagnostic use (IVD). Please refer to manufacturer's product insert for further information. For product use outside of that specified on the manufacturer's labelling and/or instructions for use, Bio-Strategy does not take responsibility.