

-30°C Chest Freezers

PHCbi

As a storage environment with excellent safety features, easy operability, and a host of other features, these freezers offer unsurpassed reliability and functionality. If you are looking for precision temperature-controlled storage equipment as low as -30°C, trust the PHCbi Biomedical Freezer for your valuable samples.

Only \$2,737

Usually \$3,910

Code: SANYMDF-137-PE
138L with 2 baskets



**Limited Stock
Order Today!**

Only \$3,079

Usually \$4,399

Code: SANYMDF-237-PE
221L with 3 baskets



- Temperature range: -20°C to -30°C
- Hi/Lo temperature & power failure alarms & dry remote alarm contacts
- Quiet in operation: 41dB (A)
- Hermetic compressor with low energy consumption
- Alarms for high and low temperature, door open, and power failure
- Two access ports
- User-friendly controller with digital display.

- Display keylock prevents adjustment of the control panel by unauthorised personnel
- Compressor delay timer for staggered start-up after power outage, preventing circuit overload
- Castors for mobility and locating feet when in working position
- Keyed lid lock for security
- Suitable for operation on 220/240V 50Hz 10 amp supply
- 5-year warranty

Valid Until: 28 Feb 2026

Ref. Code: SANY-Chest

0800 34 24 66 | sales.nz@bio-strategy.dksh.com | www.bio-strategy.com | labshop.dksh.co.nz

Prices and discounts may be subject to change and only available while stocks last. Not valid with any other special offer or discounts and cannot be applied in retrospect. All prices in local currency excluding GST and local delivery charges. Products advertised may be for research use only (RUO) or in vitro diagnostic use (IVD). Please refer to manufacturer's product insert for further information. For product use outside of that specified on the manufacturer's labelling and/or instructions for use, Bio-Strategy does not take responsibility.

Delivering Growth

 **bio-strategy**
Part of DKSH Group

 **DKSH**