

# Gilson Product Brochure

Gilson provides high quality, dependable solutions for today's demanding liquid chromatography, solid phase extraction, liquid handling and gel permeation chromatography clean-up requirements. From the legendary PIPETMAN® products to a full range of system solutions, Gilson is the element of purification for your laboratory.

## Air Displacement Pipettes



### PIPETMAN/PIPETMAN L

Recognized as the pipetting standard for its reliability and proven durability, PIPETMAN® combines legendary accuracy, precision, and durability with easy pipetting and handling.

### MyPIPETMAN Select

Ergonomic, safe, and reliable, MyPIPETMAN® Select offers volume setting choice and optional handle personalization for optimal use in various protocols.

## Pipette Tips



### Gilson Pipette Tips

Designed to fit perfectly and provide the best results with Gilson pipettes, the universal design of Gilson pipette tips allows them to be used with most leading pipette brands.

## Positive Displacement Pipettes



### MICROMAN E

MICROMAN® E positive-displacement pipettes, along with capillary pistons, offer the highest precision when ergonomically pipetting non-aqueous samples, such as viscous or volatile liquids.

### DISTRIMAN

DISTRIMAN® is a positive-displacement, continuously adjustable, repetitive pipette with direct read out designed to simplify multiple dispensing when used with DISTRITIPS® syringes.

## Electronic Pipettes



### PIPETMAN M Connected

PIPETMAN® M Connected is a Bluetooth®-enabled electronic pipette that offers accuracy and precision in both standard and repetitive pipetting modes.

## Liquid Handling



### MACROMAN

MACROMAN® is a reliable pipette controller for glass and plastic serological pipettes from 1 to 100 mL.

MACROMAN aspiration and dispensing speed can be set using pressure sensitive triggers or by using the thumbwheel while pipetting. The speed is adjustable up to 25 mL in less than three seconds. The screen continuously displays battery status and the aspiration and dispensing speed setting.

## Centrifuges



### CENTRY 103 Minicentrifuge

CENTRY™ 103 Minicentrifuge is an easy-to-use and portable minicentrifuge, ideal for the quick separation of particles from a supernatant or for driving down of a drop of liquid from the walls of microtubes. Thanks to its ergonomic design, small footprint, and optional battery control, the CENTRY 103 can be used anywhere.

## Benchtop Instruments



### Safe Aspiration Station

The Gilson Safe Aspiration Station is a stand-alone motorized vacuum system for quick and easy aspiration of biological liquid waste disposal. Liquids are quickly and safely aspirated and collected in a 4 L polypropylene bottle. The station offers adjustable aspiration speed and handles, ideal for ELISA tests, cell culture, and any other contaminated or hazardous liquid aspiration and disposal.

## Benchtop Instruments



### DISPENSMAN

DISPENSMAN® is a bottle-top dispenser for convenient, safe, and precise delivery of common laboratory solutions and solvents. Its innovative, three-position nozzle controls the dispensing of liquids, the removal of air in the system to avoid liquid loss, and includes an anti-drip safety feature that returns liquid to the bottle. Make volume changes quickly and easily with no need to open bottles, decreasing chances of contamination.



### Mini Vortex Mixer

The Mini Vortex Mixer delivers a fixed speed of 2,800 rpm by only applying slight pressure to the center of the spill-resistant vortex head. Remove the pressure and the vortex motor immediately stops. The entire operation is quiet and smooth and contained in a small 10 x 10 cm (4 x 4 in.) footprint.



### Roto-Mini Plus Rotator

The Roto Mini Plus variable speed rotator offers three operational mixing modes and accommodates many standard tube sizes, making it a flexible addition to any lab.



labshop.dksh.co.nz

## DKSH LabShop

Register an account with us today to start shopping lab products online.

[labshop.dksh.co.nz](http://labshop.dksh.co.nz)

## Contact Us

Phone: 0800 951 010

Email: [customerservicenz@dksh.com](mailto:customerservicenz@dksh.com)

Address: 94 Apollo Drive, Rosedale, Auckland 0632, New Zealand