

Axygen® PCR Microplates Compatibility Chart

		STANDARD PCR THERMAL CYCLERS																																																
		Axygen®	Applied Biosystems®					Eppendorf®											BioRad®/MJ Research®							Agilent®	Analytik Jena/Biometra																							
Corning Cat. No.		MaxyGene II	2720	9600	9700	9800 Fast	Veriti®/Fast	Mastercycler®	Mastercycler ep	Mastercycler M384	Mastercycler nexus	Mastercycler nexus eco	Mastercycler nexus flat	Mastercycler nexus flat eco	Mastercycler nexus gradient	Mastercycler nexus gradient eco	Mastercycler nexus GSX1	Mastercycler nexus GSX1e	Mastercycler nexus SX1	Mastercycler nexus SX1e	Mastercycler pro 384	Mastercycler pro/S	Cl,000™ Touch™	DNA Engine Dyad®/Dyad Disciple™	Engine Tetrad® 2	MyCycler™	PTC-100®	PTC-200	PTC-225 Tetrad	S1,000™	T100™	SureCycler® 8800	FlexCycler	T1 Thermal Cycler	TGradient	TPersonal	TProfessional/Basic	TProfessional/Standard	TProfessional/TRIO	TRobot	Uno	Uno II								
PCR Microplate	PCR-96-FLT-C	■	■	■	■	■	■	■	■	■	■	■	■	■	■	■	■	■	■	■	■	■	■	■	■	■	■	■	■	■	■	■	■	■	■	■	■	■	■	■	■	■	■	■	■	■	■			
	PCR-96-LP-FLT-C	■				■	■*																																											
	PCR-96-C		■	■	■		■								■	■	■	■	■	■	■	■	■																											
	PCR-96-MB-C																																																	
	PCR-96-SG-C		■	■	■		■																	■	■	■	■	■	■	■	■	■	■	■	■	■	■	■	■	■	■	■	■	■	■	■	■	■		
	PCR-96M2-HS-C	■	■	■	■		■	■	■			■	■			■	■	■	■	■	■	■	■	■	■	■	■	■	■	■	■	■	■	■	■	■	■	■	■	■	■	■	■	■	■	■	■	■	■	
	PCR-96-HS-AC-C	■	■	■	■		■	■	■			■	■			■	■	■	■	■	■	■	■	■	■	■	■	■	■	■	■	■	■	■	■	■	■	■	■	■	■	■	■	■	■	■	■	■	■	
	PCR-96-HS-C	■	■	■	■		■	■	■			■	■			■	■	■	■	■	■	■	■	■	■	■	■	■	■	■	■	■	■	■	■	■	■	■	■	■	■	■	■	■	■	■	■	■	■	
	PCR-96-AB-C		■	■	■		■																																											
	PCR-96-LP-AB-C					■	■*																																											
	PCR-96-FS-C	■						■	■			■	■			■	■	■	■	■	■	■	■	■	■	■	■	■	■	■	■	■	■	■	■	■	■	■	■	■	■	■	■	■	■	■	■	■	■	■
	PCR-96-LC480-W-NF																																																	
	PCR-384-C				■		■			■												■	■	■	■	■	■	■	■	■	■	■	■	■	■	■	■	■	■	■	■	■	■	■	■	■	■	■	■	
	PCR-384M2-C				■		■			■											■	■	■	■	■	■	■	■	■	■	■	■	■	■	■	■	■	■	■	■	■	■	■	■	■	■	■	■	■	
	PCR-384-RGD-C				■		■			■											■	■	■	■	■	■	■	■	■	■	■	■	■	■	■	■	■	■	■	■	■	■	■	■	■	■	■	■	■	
	PCR-384-LC480-W-NF																																																	
PCR Strip/ Individual Tube	PCR-0208-C	■	■	■	■		■	■			■	■			■	■	■	■	■	■	■	■	■	■	■	■	■	■	■	■	■	■	■	■	■	■	■	■	■	■	■	■	■	■	■	■	■	■	■	■
	PCR-0108-LP-C	■				■	■*																																											
	PCR-0108-LP-RT-C					■	■*																																											
	PCR-0212-C	■	■	■	■		■	■				■	■			■	■	■	■	■	■	■	■	■	■	■	■	■	■	■	■	■	■	■	■	■	■	■	■	■	■	■	■	■	■	■	■	■	■	■
	PCR-02-C	■	■	■	■		■	■				■	■	■	■	■	■	■	■	■	■	■	■	■	■	■	■	■	■	■	■	■	■	■	■	■	■	■	■	■	■	■	■	■	■	■	■	■	■	

* Compatible with ABI® Fast PCR only.

