

Product information



CULTURA® Mini-Incubator

Art. No. ZC7140651

The CULTURA® Mini-Incubator has been used for 20 years in the fields of medical and dental diagnostics, food bacteriology and hygiene monitoring, monitoring of coolants, lubricants and industrial liquids as well as of drinking water, industrial process and swimming pool water. The handy small incubator is suitable for the incubation of RIDA®STAMP, Compact Dry, VitaFast® microtiter plates as well as for all RIDASCREEN® ELISA, where an incubation step at 35 - 37 °C is intended. Particularly small and middle sized premises can use it to install a cost-effective lab.



CULTURA® Mini-Incubator, Art Nr. ZC7140651

The incubator is suitable for all incubations between 25 to 45 °C. The thermostat can be adjusted via a setscrew at the front side of the cabinet utilizing the thermometer to control the temperature setting. Further advantages of the device are its compact, stable housing, easy cleaning and the full-view plexiglass door, which allows visual inspection at any time.

Technical Data			
Trade name:	CULTURA®	Door	PC transparent
Art. No.	ZC7140651	Thermometer	Up to 60 °C, in special tube
Temperature range	25 - 45 °C	Base cover	ABS
Temperature stability	+/- 1 °C	Warranty	2 years from initial installation
Voltage	230 V/110 V/12 V	External dimensions (W x H x T)	310 x 155 x 168 mm
Frequency	50 Hz	Internal dimensions (W x H x T)	220 x 120 x 150 mm
Power output	26 W/0.2 A/220 V	Weight	1.1 kg
Protective class	II		
Case	Not easily flammable PC & ABS		



bio-strategy
delivering technology

in **New Zealand** please contact

Bio-Strategy Ltd

T: 0800 34 24 66 | E: sales.nz@bio-strategy.com

www.bio-strategy.com | shop.bio-strategy.com