

Scharlab Reader App for iOS and Android



Scharlab Reader App lets you obtain technical documents and mini safety labels by simply scanning the QR code on your Scharlau bottle.

Instantly, at any time, users can download mini safety labels, Certificates of Analysis (CoA), Technical Data Sheets (TDS) and Safety Data Sheets (SDS) for their Scharlau products.

- The SDS are available in 16 languages (English, Spanish, German, French, Italian, Chinese, Polish, Danish, Portuguese, Australian English, Finnish, Greek, Hungarian, Norwegian, Romanian, Swedish) and the mini safety labels in 9 languages (English, Spanish, German, French, Italian, Chinese, Polish, Danish, Portuguese).
- The mini safety labels include colour hazard symbols, hazard indications (H) and precautionary statements (P). Since they are downloaded from our website, the user can always access the most updated versions.
- The user interface is available in 6 languages (Spanish, English, Italian, French, German and Portuguese).



Mini safety label

Art. No.

Description

Hazard statements (H)

Precautionary statements (P)

Colour hazard symbols

MEQ3151000


Methanol, Multisolvant® HPLC grade ACS ISO UV-VIS K.F.

Highly flammable liquid and vapour. Toxic if swallowed. Toxic in contact with skin. Toxic if inhaled. Causes damage to organs.

Keep away from heat / sparks / open flames / hot surfaces. - No smoking. Use explosion-proof electrical / ventilating / lighting / equipment. IF ON SKIN (or hair): Remove / Take off immediately all contaminated clothing. Rinse skin with water / shower. Remove / Take off immediately all contaminated clothing. Store locked up. Dispose of contents / container in accordance with local / regional / national / international regulations.

You can find the application on Apple Store and Google Play and install it like any other App.

1. The App welcome screen



2. The App main screen.

With the "Scan" button you access the screen where all the readings are collected.




3. Make a new reading

Press the [+] button to start the camera to make a new reading.




To cancel the scan, select the grey square that is below the image.

4. If you want to delete a previous reading, slide your finger to the left on the reading you wish to eliminate. A "Delete" button will appear, which will remove the item.




5. When pressing on a reading, it leads you to a screen where you can access several available documents.




6. Selecting the type of document

Pressing the dropdown menu button, you can select among several documents (labels, SDS, TDS, COA).




7. Printing the documents

After pressing "Download", the selected document will start downloading. If the download is correct, an icon will appear.




8. Printing the documents

After pressing on the document icon, it will be shown on the screen.




9. Sending the documents

After pressing the "Send to" button, the mail application will open. The user can easily send the document to any email.



10. Introducing email address

Click on the "Settings" button to introduce the email address where you wish to send the documents. You can also choose the language of the App in the settings screen.



11. "Info"

The button "Info" shows the terms and conditions of the App.

