

15% off Axygen Equipment



Valid Until:
31 December 2018

Ref. Code:
AXYG-EQUIP

15%
OFF

Quality lab equipment from Axygen now on sale - talk to your Sales Specialist today!

Speak to your local Bio-Strategy Representative for further details or call Customer Service on 0800 34 24 66 or email sales.nz@bio-strategy.com

Axygen® Axyspin Mini Plate Spinner Centrifuge

The Axygen Axyspin Plate Spinner is the first centrifuge designed specifically for quick and easy centrifugation of samples in PCR plates or microplates.

- Capacity for two standard plates
- Top loading design
- Exceptionally small footprint
- Starts and stops in seconds



Purchase an Axygen Axyspin in 2018 and also receive a free timer

FREE



Axygen® Gel Documentation Systems

Axygen® Gel Documentation systems easily capture publication quality, 16 bit TIFF images. The systems are quick to set up and have an intuitive user interface for image capture, annotation, and contrast adjustment. Images are easily saved and opened in common gel analysis software for more detailed analysis.

- Auto exposure – you never have to manually manipulate the camera
- Ready-to-use – just connect the computer* and load the software
- Long-life dual wavelength transilluminators
- Exceptionally high resolution (5.4 MP) camera
- The cabinet is compact with an ultra-violet (UV) cut-off switch to prevent exposure to UV light if the cabinet door is opened while the UV light is on



Axygen® MaxyGene II Thermal Cycler

The Axygen MaxyGene II Thermal Cycler brings increased speed and advanced features, providing the premium performance you have come to expect from Axygen brand products.

- Unique flexible programming
- Rapid run times
- Improved workflow over traditional gradient cyclers
- Ramping rates up to 5°C/sec
- Simple user interface
- Adjustable heated lid accommodates strips, tubes, and microplates
- Pre-programmed standard protocols



More equipment on the next page...

BIO-STRATEGY LIMITED

Tel: 0800 34 24 66
Email: sales.nz@bio-strategy.com
www.bio-strategy.com



bio-strategy
delivering technology

Prices & discounts may be subject to change and only available while stocks last. Not valid with any other special offer or discounts and cannot be applied in retrospect. All prices exclude shipping & delivery charges. Promotion code (if stated) must be referenced on your purchase order to receive offer prices.

Axygen[®] Horizontal Gel Boxes

Axygen horizontal gel boxes are high quality, full featured units designed for exceptional durability and efficiency. The gel boxes feature leak-free external gel casting, which allows you to cast one gel outside of the box while another gel is running inside.



- Flexibility – each unit comes with 2 sizes of gel trays to optimise applications
- External gel casting to gel box allows for running gels while casting the next one
- The Axxygen contrast enhancer makes loading gels easier
- Rugged injection molding is supported by a 3-year warranty
Larger units come with recirculation ports allowing for buffer exchange

Axygen[®] PlateMax[®] Semi-automatic Plate Sealer

The Axxygen PlateMax semi-automatic plate sealer is ideal for the low to medium throughput laboratory that requires uniform and consistent sealing of microplates. Sealing parameters are set and displayed via the user-friendly control panel and the sealing operation is automated to ensure consistent results.



- Accepts a wide range of microplates and heat seals
- User-friendly control system
- Large digital display shows all operating parameters
- Adjustable sealing temperature
- Precise temperature, timing, and pressure for consistent sealing

Axygen[®] Axyspin Refrigerated Microcentrifuge

The Axxygen Axyspin refrigerated microcentrifuge is a high performance refrigerated microcentrifuge designed to quickly process samples for applications such as nucleic acid and protein preparation, purifications, extractions, and pelleting.



- Powerful refrigeration system cools to 4°C in 8 minutes
- 24 x 1.5/2.0 mL aluminium rotor included
- Exceptionally quiet and compact

Talk to your Sales Specialist today!

**15%
OFF**

BIO-STRATEGY LIMITED

Tel: 0800 34 24 66
Email: sales.nz@bio-strategy.com
www.bio-strategy.com



bio-strategy
delivering technology