

Valid Until:  
31 December 2018

Ref. Code:  
CORN-COOL

## Savings on CoolCell LX range

CoolCell LX, in combination with a -80°C freezer, will provide the freezing rate of -1°C/minute that is ideal for cryopreservation of most cells and cell lines. Using a combination of closed cell cross-linked foam housing, a thermo-conductive alloy core, and radial vial symmetry, freezing profiles are consistent and reproducible. The low heat content also ensures that CoolCell LX containers will rapidly return to room temperature when removed from the freezer.

Unique features of controlled-rate freezing with Corning CoolCell® include:

- Ease of use
- Alcohol and fluid-free freezing
- Lower cost of use than alcohol-based devices
- High cell recovery and cell viability
- Reproducibility
- Simple, consistent way to standardise controlled-rate freezing

The Corning Cryogenic Storage Solutions ensure standardised, controlled-rate freezing for all cell types. The Corning sample cooling and heating systems ensures temperature consistency and minimize contamination risk.

**Speak to your local Bio-Strategy Representative for further details or call Customer Service on 0800 34 24 66 or email [sales.nz@bio-strategy.com](mailto:sales.nz@bio-strategy.com)**



Description	Pack Size	Code	Sale Price (excl. GST)
Corning® CoolCell® LX, Cell Freezing Container, for 12x1 mL or 2mL Cryogenic Vials, Purple	1 item	CORN432001	<b>\$248.00</b>
Corning® CoolCell® LX, Cell Freezing Container, for 12x1 mL or 2mL Cryogenic Vials, Green	1 item	CORN432002	<b>\$248.00</b>
Corning® CoolCell® LX, Cell Freezing Container, for 12x1 mL or 2mL Cryogenic Vials, Orange	1 item	CORN432003	<b>\$248.00</b>
Corning® CoolCell® LX, Cell Freezing Container, for 12x1 mL or 2mL Cryogenic Vials, Pink	1 item	CORN432004	<b>\$248.00</b>

BIO-STRATEGY LIMITED

Tel: 0800 34 24 66  
Email: [sales.nz@bio-strategy.com](mailto:sales.nz@bio-strategy.com)  
[www.bio-strategy.com](http://www.bio-strategy.com)



Prices & discounts may be subject to change and only available while stocks last. Not valid with any other special offer or discounts and cannot be applied in retrospect. All prices exclude shipping & delivery charges. Promotion code (if stated) must be referenced on your purchase order to receive offer prices.