

# DURAN® UTILITY Laboratory Bottles



Valid Until:  
31 December 2018

Ref. Code:  
SCOT-YOU

## DURAN® YOUTILITY clear and amber Laboratory bottles with GL 45 thread

Clear and amber DURAN® YOUTILITY bottles with slimmer shape for optimal space usage in your autoclaves and refrigerators.

- Manufactured from Duran® borosilicate 3.3 pharmacopoeia type 1 neutral glass
- Ergonomic designed cap for easier tightening or removing especially for gloved or small hands
- New GL45 din thread for faster and easier opening
- GL45 coloured tags for easy identification of bottles in your laboratory
- Durable white ceramic markings for clear and accurate reading of the volumes
- Specially shaped gripping zones on each side for easier and safer handling
- Amber bottles for storage or transport of light sensitive chemicals from uv light
- Amber bottles meet the requirements of usp <660> and eur. ph. 3.2.1

**Speak to your local Bio-Strategy Representative for further details or call Customer Service on 0800 34 24 66 or email [sales.nz@bio-strategy.com](mailto:sales.nz@bio-strategy.com)**



Description	Unit of Measure	Code	List Price (excl. GST)	Sale Price (excl. GST)
Bottle 125ml Youtility clear w/cyan GL45 screw cap	4 items	SCOT218812854	\$45.83	\$34.37
Bottle 250ml Youtility clear w/cyan GL45 screw cap	4 items	SCOT218813653	\$51.74	\$38.80
Bottle 500ml Youtility clear w/cyan GL45 screw cap	4 items	SCOT218814452	\$57.88	\$43.41
Bottle 1L Youtility clear w/cyan GL45 screw cap	4 items	SCOT218815457	\$80.79	\$60.59
Bottle 125ml Youtility amber w/cyan GL45 screw cap	4 items	SCOT218862859	\$204.90	\$153.67
Bottle 250ml Youtility amber w/cyan GL45 screw cap	4 items	SCOT218863658	\$244.87	\$183.65
Bottle 500ml Youtility amber w/cyan GL45 screw cap	4 items	SCOT218864457	\$259.13	\$194.34
Bottle 1L Youtility amber w/cyan GL45 screw cap	4 items	SCOT218865453	\$366.11	\$274.58

25%  
OFF



BIO-STRATEGY LIMITED

Tel: 0800 34 24 66  
Email: [sales.nz@bio-strategy.com](mailto:sales.nz@bio-strategy.com)  
[www.bio-strategy.com](http://www.bio-strategy.com)



Prices & discounts may be subject to change and only available while stocks last. Not valid with any other special offer or discounts and cannot be applied in retrospect. All prices exclude shipping & delivery charges. Promotion code (if stated) must be referenced on your purchase order to receive offer prices.