

Falcon® Culture Slides

FALCON®

A Corning Brand

Valid Until:
31 December 2018

Ref. Code:
BDAA-SLIDE

Innovative design, safe and easy to use!

Falcon Culture Slides allow you to culture cells and then analyse them on a glass microscope slide. Cells are grown in a plastic chamber affixed to a specially prepared glass microscope slide. Cells can be fixed and stained in place without disruption of the cell monolayer. The chamber is easily and safely removed with an easy-to-use, disposable Safety Removal Tool.

Designed for consistent cell culture results

- Specially cleaned and triple-rinsed glass slides
- Performance validated with HEP-2 and BAE cells
- Tested for 72-hour confluency with MRC-5 and BAE cells

Features

- 1.2mm bevelled-edge slide, 25 mm x 75 mm, soda-lime glass
- Pressure-sensitive, biocompatible, non-migrating, acrylic-adhesive gasket
- Blue hydrophobic border defines cell culture areas
- Polystyrene vessel, lid, and tool
- Supplied with disposable Safety Removal Tool
- Wells numbered for easy identification
- Sterile
- Trays designed for use in incubator

Speak to your local Bio-Strategy Representative for further details or call Customer Service on 0800 34 24 66 or email sales.nz@bio-strategy.com

Pack Size	Code	Sale Price (excl. GST)
Falcon® 4-well culture slide, glass slide with polystyrene vessel, lid, and safety removal tool		
96 items	BDAA354104	\$856.80
96 items	BDAA354108	\$856.80
24 items	BDAA354114	\$267.60
24 items	BDAA354118	\$278.46

40%
OFF!



BIO-STRATEGY LIMITED

Tel: 0800 34 24 66
Email: sales.nz@bio-strategy.com
www.bio-strategy.com



bio-strategy
delivering technology

Prices & discounts may be subject to change and only available while stocks last. Not valid with any other special offer or discounts and cannot be applied in retrospect. All prices exclude shipping & delivery charges. Promotion code (if stated) must be referenced on your purchase order to receive offer prices.