

Integra Adjustable Spacing Pipettes Promotion

INTEGRA

Valid Until:
31 December 2018

Ref. Code:
INTG-TRADE

Trade-in any adjustable tip spacing pipette and save on a new VOYAGER!

Trade-in any multichannel pipette and save on a new VOYAGER!

With the push of a button, perform parallel transfer of multiple samples between labware of different sizes and formats.


VOYAGER multichannel pipettes have been designed to address the problem of single transfers such as gel loading, tube to plate transfers or 96-well to 384-well plate reformatting.

It is the first and only automatic adjustable tip spacing pipette on the market. Using a VOYAGER you can transfer multiple samples at once in the blink of an eye. Moreover, by reducing transfer steps, VOYAGER helps you prevent pipetting errors and repetitive strain injuries (RSI).

The motorised tip spacing offers many advantages over manual tip spacing adjustments:

- Single handed operation leaves the other hand free to handle labware
- On the fly access to up to three user-defined tip spacings
- These user defined tip spacings are saved, no need to memorise spacings of different labware formats

Speak to your local Bio-Strategy Representative for further details or call Customer Service on 0800 34 24 66 or email sales.nz@bio-strategy.com



BIO-STRATEGY LIMITED

Tel: 0800 34 24 66
Email: sales.nz@bio-strategy.com
www.bio-strategy.com

 **bio-strategy**
delivering technology

Prices & discounts may be subject to change and only available while stocks last. Not valid with any other special offer or discounts and cannot be applied in retrospect. All prices exclude shipping & delivery charges. Promotion code (if stated) must be referenced on your purchase order to receive offer prices.