

HistoCore MULTICUT Semi-motorised Rotary Microtome

Leica
BIO SYSTEMS

Valid Until:
31 December 2018

Ref. Code:
LEIB-HISTO

The NEXT generation with Leica Proven technology

What's new?

Personal coarse feed wheel

- Manual or Motorised use with the push of a button

Fast home function

- 2x faster than previous model (RM2245)

Programmable MEMO position

- Rapid specimen exchange while trimming
- Remembers the optimal specimen position for your sectioning

Antistatic Waste tray

- Shorter cleaning times from minutes to seconds

40% more flat space on microtome top

- Allows placement of larger objects on top

Speak to your local Bio-Strategy Representative for further details or call Customer Service on 0800 34 24 66 or email sales.nz@bio-strategy.com

NEW!



BIO-STRATEGY LIMITED

Tel: 0800 34 24 66
Email: sales.nz@bio-strategy.com
www.bio-strategy.com



bio-strategy
delivering technology

Prices & discounts may be subject to change and only available while stocks last. Not valid with any other special offer or discounts and cannot be applied in retrospect. All prices exclude shipping & delivery charges. Promotion code (if stated) must be referenced on your purchase order to receive offer prices.