

# Surgipath® Paraplast® Tissue Embedding Medium



Valid Until:  
31 December 2018

Ref. Code:  
LEIS-PARA

## Different pack size = lower price per kg

**Paraplast is recommended for general tissue embedding.**

- Melting point of 56°C
- Effective tissue infiltration
- Superior quality sections with minimal compression
- Gives brittle, difficult to section tissue an innate elasticity, promoting wrinkle-free sections

**Speak to your local Bio-Strategy Representative for further details or call Customer Service on 0800 34 24 66 or email [sales.nz@bio-strategy.com](mailto:sales.nz@bio-strategy.com)**

Description	Pack Size	Code	List Price (excl. GST)
Paraplast Regular 5kg bags (x4)	4x 5kg bags (20kg)	LEIS39602012	<b>\$189.49</b>

Price per kg: \$9.47/kg.



**Ex-stock**

BIO-STRATEGY LIMITED

Tel: 0800 34 24 66  
Email: [sales.nz@bio-strategy.com](mailto:sales.nz@bio-strategy.com)  
[www.bio-strategy.com](http://www.bio-strategy.com)



Prices & discounts may be subject to change and only available while stocks last. Not valid with any other special offer or discounts and cannot be applied in retrospect. All prices exclude shipping & delivery charges. Promotion code (if stated) must be referenced on your purchase order to receive offer prices.