

Marienfeld Slides and Coverslips Now on Sale!



Valid Until:
31 December 2018

Ref. Code:
MARI-COVER

**SPECIAL
PRICES**

SPECIAL PROMO PRICES VALID UNTIL END OF 2018

The brand Marienfeld Superior is a synonym for quality and reliability, with more than 90 years experience in developing and producing laboratory glassware



Marienfeld Superior microscope slides stand out by their exceptional wettability and cleanliness. They are manufactured to exact size and are therefore well suited for application in automatic machinery.

- Made of soda lime glass of third hydrolytic class
- In compliance with DIN ISO 8037/1
- Thickness: approx. 1mm (tol. $\pm 0.05\text{mm}$) - other thicknesses available
- Pre-cleaned and ready to use
- Frosted microscope slides have a silky, frosted marking area of approx. 20mm on both sides

Marienfeld Superior cover glasses feature an excellent flatness and smooth surface thanks to the best material. No cover glasses made of soda lime glass match the quality of Marienfeld borosilicate glass coverslips.

- Absolutely colourless, perfectly clear, suitable for fluorescence microscopy
- In compliance with DIN ISO 8255

Speak to your local Bio-Strategy Representative for further details or call Customer Service on 0800 34 24 66 or email sales.nz@bio-strategy.com



Description	Unit of Measure	Code	Sale Price (excl. GST)
Slide 75x26mm Twin Frost Cut Edge	50	MARI1005200	\$3.95
Slide 75x26mm Plain Ground Edge	50	MARI1005412	\$4.95
Slide 75x26mm Twin Frost Ground Edge	50	MARI1005612	\$5.65
Coverslips 22x22mm No.1	100	MARIO101050	\$3.75
Coverslips 24x40mm No 1	100	MARIO101192	\$3.75
Coverslips 18x18mm No.1	100	MARIO101030	\$2.45
Coverslips 22x40mm No.1	100	MARIO101122	\$6.50
Coverslips 22x50mm No.1	100	MARIO101142	\$6.95
Coverslips 22x60mm No.1	100	MARIO101152	\$7.30
Coverslips 24x50mm No.1	100	MARIO101222	\$7.50
Coverslips 24x60mm No.1	100	MARIO101242	\$8.90
Coverslips 15 x15mm No.0	100	MARIO100010	\$1.95

For our Leica Microscope users we are offering a free pack at no charge (maximum 1 per microscope, maximum 5 per lab). Contact us to claim your free pack today!

BIO-STRATEGY LIMITED

Tel: 0800 34 24 66
Email: sales.nz@bio-strategy.com
www.bio-strategy.com



Prices & discounts may be subject to change and only available while stocks last. Not valid with any other special offer or discounts and cannot be applied in retrospect. All prices exclude shipping & delivery charges. Promotion code (if stated) must be referenced on your purchase order to receive offer prices.