

ImageXpress Pico Automated Cell Imaging System Lab Start-Up Package

Valid Until:
14 December 2018

Ref. Code:
MDEV-IXPLS



Purchase before 14 December and receive a lab start-up configuration for under \$100,000

Intelligent cell imaging and analysis integrated in a single, affordable system

More than a digital microscope, the ImageXpress® Pico system combines imaging and powerful analysis for individual labs that need easy automated imaging solutions. With over 25 preconfigured cell imaging protocols available, the ImageXpress Pico system can be used for a variety of cell biology assays. With a simple step-by-step setup, images are automatically acquired and analysed while users can work on other tasks.

Get imaging quickly with the Lab Start-Up Package:


- Leica FLUOTAR 4x/NA 0.13 Objective
- Leica FLUOTAR 10x/NA 0.32 Objective
- Transmitted light capability
- FITC filter cube for fluorescence imaging
- Microtiter plate adapter and 4-slide adapter
- CellReporterXpress Software license – WITH 5 Logins from a desktop AND 17 analysis modules (cell counting, scoring, apoptosis to name a few)
- System includes 1-year warranty covering parts & labour

Speak to your local Bio-Strategy Representative for further details or call Customer Service on 0800 34 24 66 or email sales.nz@bio-strategy.com

**additional objectives, filters and modules are available for purchase.*



**Save over
\$20,000**



BIO-STRATEGY LIMITED

Tel: 0800 34 24 66
Email: sales.nz@bio-strategy.com
www.bio-strategy.com



Prices & discounts may be subject to change and only available while stocks last. Not valid with any other special offer or discounts and cannot be applied in retrospect. All prices exclude shipping & delivery charges. Promotion code (if stated) must be referenced on your purchase order to receive offer prices.