

PHCbi MDF-U731M -30°C Upright Biomedical Freezer, 690 Litres



Valid Until:
31 December 2018

Ref. Code:
MDF-U731M

Direct cooling and manual defrost for storing temperature sensitive samples

Large capacity with 4 adjustable wire shelves

The PHCbi MDF-U731M Upright Biomedical -30°C Freezer has a total capacity of 690 litres. The unit has manual defrost and is fully alarmed (high/low temperature, power failure, door ajar & compressor temperature alarms). The 4 adjustable wire shelves can support up to 50kg/shelf. The unit is very quiet in operation, runs off a standard 10 amp supply.

- Powerful refrigeration system achieves -30°C at up to a 35°C ambient
- User friendly eye-level digital display with audible and visual alarms
- Large-capacity, economical and space-saving model
- Very quiet in operation: 42dB(A)
- Keyed door lock for security and manual defrost to prevent temperature spiking
- Dry remote alarm contacts as standard

Speak to your local Bio-Strategy Representative for further details or call Customer Service on 0800 34 24 66 or email sales.nz@bio-strategy.com

Description	Code	List Price (excl. GST)	Sale Price (excl. GST)
-30°C Upright Biomedical Freezer, 690 litres	SANYMDF-U731M	\$8,991.00	\$7,995.00



SALE

BIO-STRATEGY LIMITED

Tel: 0800 34 24 66
Email: sales.nz@bio-strategy.com
www.bio-strategy.com



bio-strategy
delivering technology

Prices & discounts may be subject to change and only available while stocks last. Not valid with any other special offer or discounts and cannot be applied in retrospect. All prices exclude shipping & delivery charges. Promotion code (if stated) must be referenced on your purchase order to receive offer prices.