

# PHCbi MPR-514 Upright Pharmaceutical Refrigerator, 2°C to 14°C, 489 Litres

Valid Until:  
31 December 2018

Ref. Code:  
MPR-514



PHCbi High Performance  
Sliding Door Pharmaceutical  
Refrigerators.

## Smaller footprint due to dual sliding doors – but capacity of 489 litres!

The PHCbi MPR-514 series of Upright Pharmaceutical Refrigerators are available in two models; MPR-514 with 5 plastic coated wire shelves (50kg load/shelf) or the MPR-514R with 5 plastic coated wire shelves (20kg load/shelf) and 5 plastic coated sliding racks (20kg/rack).

The MPR-514 has a total capacity of 489 litres whilst the MPR-514R has a total capacity of 486 litres. Their exacting performance make the units ideal for storage of samples and reagents in pharmaceutical and research applications.

- Dual sliding door design maximises capacity with smaller operational footprint
- Fully automatic cyclic defrosting with minimal chamber temperature increase
- User friendly eye-level digital display with audible and visual alarms
- Very quiet in operation: 42dB(A)
- Keyed door lock for security
- Glass doors have protective film for safety

*In April 2018, Panasonic Healthcare Co., Ltd. changed company name to PHC Corporation. PCHbi is the new brand of PHC Corporation's biomedical division.*

**Speak to your local Bio-Strategy Representative for further details or call Customer Service on 0800 34 24 66 or email [sales.nz@bio-strategy.com](mailto:sales.nz@bio-strategy.com)**

Description	Code	List Price (excl. GST)	Sale Price (excl. GST)
Refrigerator, Pharmaceutical, 489L, 2°C to 14°C, 2 sliding doors + 5 shelves	SANYMPR-514	\$7,727.00	\$6,975.00
Refrigerator, Pharmaceutical, 486L, 2°C to 14°C, 2 sliding doors + 5 shelves + 5 racks	SANYMPR-514R	\$8,411.00	\$7,595.00

# PHCbi

ON  
SALE!



BIO-STRATEGY LIMITED

Tel: 0800 34 24 66  
Email: [sales.nz@bio-strategy.com](mailto:sales.nz@bio-strategy.com)  
[www.bio-strategy.com](http://www.bio-strategy.com)



bio-strategy  
delivering technology

Prices & discounts may be subject to change and only available while stocks last. Not valid with any other special offer or discounts and cannot be applied in retrospect. All prices exclude shipping & delivery charges. Promotion code (if stated) must be referenced on your purchase order to receive offer prices.