

Discover the Next-Generation PCR Technologies in the New QIAGEN AllTaq™

Valid Until:
31 December 2018

Ref. Code:
QIAG-TAQ

Special 40% introductory discount!

Simplify the Way You Amplify

Use the same 45-minute protocol for all your templates!

The PCR innovations in our new AllTaq have eliminated annealing temperature and buffer composition optimisation and allows confirmation of correct pipetting with a simple colour change. Reliably amplify GC-rich or long targets, up to 9kb, with sensitivity to a single target molecule.

The AllTaq Master Mix and AllTaq PCR Core Kits can be used for any PCR, for genotyping and genetic analysis or target detection from complex gDNA or cDNA templates, including reactions with multiple primer pairs or that include a wide range of amplicon sizes.

Try AllTaq and have complete confidence in your PCR!

Speak to your local Bio-Strategy Representative for further details or call Customer Service on 0800 34 24 66 or email sales.nz@bio-strategy.com

Description	Code	List Price (excl. GST)	Sale Price (excl. GST)
AllTaq PCR Core Kit (250 U)	QIAG203123	\$291.00	\$174.60
AllTaq PCR Core Kit (1000 U)	QIAG203125	\$1,014.00	\$608.40
AllTaq PCR Core Kit (5000 U)	QIAG203127	\$3,412.00	\$2,047.20
AllTaq Master Mix Kit (500)	QIAG203144	\$707.00	\$424.20
AllTaq Master Mix Kit (2500)	QIAG203146	\$2,941.00	\$1,764.60



40%
OFF!



BIO-STRATEGY LIMITED

Tel: 0800 34 24 66
Email: sales.nz@bio-strategy.com
www.bio-strategy.com



bio-strategy
delivering technology

Prices & discounts may be subject to change and only available while stocks last. Not valid with any other special offer or discounts and cannot be applied in retrospect. All prices exclude shipping & delivery charges. Promotion code (if stated) must be referenced on your purchase order to receive offer prices.